Abstract ID: 121
Presenter ID: D22.4




# IDENTIFICATION OF A LIPOLYTIC Trichoderma sp. AND CHARACTERIZATION OF ITS EXTRACELLULAR LIPASE

RESCON 2021 – ORAL PRESENTATIONS

Department of Molecular Biology & Biotechnology Faculty of Science University of Peradeniya

N.K. Athukorala, Prof. P. Samaraweera



#### **CONTENTS**

01

#### INTRODUCTION

What are lipases?
Significance of fungal lipases
Objectives

02

#### **METHODOLOGY**

Establishment of pure cultures
Screening lipolytic activity
Identification of the fungus
Optimization of growth medium
Characterization of lipase
activity

03

#### **RESULTS & DISCUSSION**

Pure cultures and lipolytic activity
Identification of the fungus
Optimized growth conditions
Optimized lipase activity

04


**CONCLUSIONS** 

### 01

**INTRODUCTION** 

#### What are lipases?

- Lipases are a significant group of biocatalysts, hydrolyzing carboxylic ester bonds to release carboxylic acids and alcohols (Daiha et al., 2015)
- Lipases catalyze other industrially demanding reactions such as
  - ✓ Esterification
  - ✓ Transesterification
  - ✓ Acidolysis (Mehta et al., 2017)
- Most of the lipases get activated by a water-lipid interface



**Figure 1:** Diagrammatic representation of the overall mechanism of lipase catalysis

#### **Significance of fungal lipases in industries**

- ✓ Low cost of production
- ✓ Availability
- ✓ Ease in genetic manipulation (Mehta *et al.*, 2017)
- ✓ Ability to tolerate polar solvents
- ✓ Thermostability
- ✓ Stability at acidic pH (Liu et al., 2015)













#### **RESEARCH OBJECTIVES**



✓ To identify the lipolytic fungus up to species-level




 To optimize the growth medium for an enhanced lipase production



To characterize lipase activity to determine optimized conditions



## **O2**METHODOLOGY



#### Identification of lipolytic activity

#### Tween 20 plate assay


#### Medium components:

Tween 20, peptone, NaCl, CaCl<sub>2</sub>, agar, and distilled water pH: 5.8

#### Phenol red plate assay

#### **Medium components:**


Phenol red indicator, olive oil, CaCl<sub>2</sub>, agar, and distilled water pH: 7.8



#### para-Nitrophenyl Palmitate (pNPP) assay (Gupta et al., 2002)

#### **Solution 1** pNPP Isopropanol Triton-X 100

**Solution 2**Gum Arabic
Tris-HCl buffer (pH 8.0)



## 03


### **RESULTS & DISCUSSION**

#### Pure cultures and lipolytic activity




Figure 2: Trichoderma colony on PDA

#### Tween 20 plate assay



**Figure 3:** White precipitate in Tween 20 plate assay

#### Phenol red plate assay




**Figure 4:** Yellow coloration in Phenol red plate assay

#### Identification of the fungus

#### Microscopic identification

- Spores: spherical (S)
- Mycelium: bidirectionally branching (M)



**Figure 5:** Microscopic image of the fungus stained with lactophenol in cotton blue


#### **Molecular identification**

- The amplified DNA fragment resulted in a band of 600 bp in 2 % agarose
- The obtained sequence showed: 100 % query cover

100 % identity


0.0 E value

with *Trichoderma longibrachiatum* sequence available in online database




**Figure 6:** The phylogenetic relationship of the unknown sequence to the selected sequences from the Gen Bank Database.


#### **Optimization of growth medium**



**Figure 7:** Affect of different carbon sources on lipase secretion from the lipolytic fungus




**Figure 8:** Affect of different nitrogen sources on lipase secretion from the lipolytic fungus




**Figure 9:** Affect of pH of the growth medium on lipase secretion from the lipolytic fungus

#### **Optimization of lipase activity**




#### Highest lipase activity at pH 6.0

Figure 10: Influenze of pH on lipase activity




#### Highest lipase activity at 40 °C

Figure 11: Influenze of temperature on lipase activity



#### Lipase activity was highest with Ca2+ as the cation

Figure 12: Influenze of different cations on lipase activity



#### Lipase activity was highest with NO<sub>3</sub>- as the anion

Figure 13: Influenze of different anions on lipase activity

# 04

**CONCLUSIONS** 

- ✓ Supported by the microscopic and molecular identification; the lipolytic fungus is *Trichoderma* longibrachiatum
- ✓ The crude lipase secretion from the lipolytic *Trichoderma* species could be increased by media optimization

Maximum lipase secretion achieved with olive oil as the carbon source, ammonium sulfate as the nitrogen source at a pH of 7.0.

✓ The crude lipase activity could be enhanced under optimized conditions

The crude enzyme activity highest at a pH of 6.0 and 40 °C. The enzyme activity could be enhanced with  $Ca^{2+}$  and  $NO_3^{-}$ .

#### **Future research**

- Further studies are required to test the enzyme's activity at high temperatures and the enzyme's thermostability
- Characterizing the enzyme activity with frequently used industrially important chemicals can expose more industrially favorable enzyme properties.



#### REFERENCES

- 1 Daiba K Angeli D. de Oliveiro C. D. & Almeide D. V. Are lineage still importe
- Daiha, K., Angeli, R., de Oliveira, S. D., & Almeida, R. V. Are lipases still important biocatalysts? A study of scientific publications and patents for technological forecasting. *PloS one*, 2015, 10(6): e0131624. DOI: https://doi.org/10.1371/journal.
- 2. Mehta, A., Bodh, U. and Gupta, R. Fungal lipases: a review. *Journal of Biotech* Research, 2017, 8(1): 58-77.
- 3. Liu, G., Hu, S., Li, L. and Hou, Y. Purification and characterization of a lipase with high thermostability and polar organic solvent-tolerance from *Aspergillus niger* AN0512. *Lipids*, 2015, **50**(11): 1155-1163. DOI 10.1007/s11745-015-4052-6.
- 4. Gupta, N., Rathi, P. and Gupta, R. Simplified para-nitrophenyl palmitate assay for lipases and esterases. *Analytical Biochemistry*, 2002, **311**(1): 98-99. DOI: 10.1016/s0003-2697(02)00379-2.

### **THANK YOU**

Any questions?

nadeeshaa@sci.pdn.ac.lk +94 77 805 1105